About Don Whitley Scientific

For over forty years we have pioneered the development and worldwide use of modified atmosphere workstations.

> More than 3000 of our workstations are in use in fifty countries.

Worldwide technical support.

Scientific guidance available from our own in-house laboratories.





Product details

The Hypoxystation is the ideal product for cell culture research. This workstation replicates conditions comparable to those occurring in vivo.

You will have the option to:

- Control O2 in 0.1% increments from 0.1% to 20% .
- Control CO2 in 0.1% increments from 0.1% to 15%
- Carry out straightforward calibration of gas sensors
- Control temperature easily and precisely
- Control relative humidity up to 80%

The Hypoxystation has been designed to combine functional reliability with operational comfort.

The touchscreen interface:

- Allows you to monitor all parameters simultaneously
- Eliminates the need for any other dials, switches and gauges

12 litre airlock provides effective cell ware transfer to and from the workstation environment in the fastest possible time

Patented, ergonomic gloveports fitted with sleeves for user comfort and ease of movement

Specification

A06070	240V	Whitley H135 HEPA Hypoxystation
A06071	110V	1452mm x 1056mm x 993m

In the interests of a policy of continuous product improvement the company reserves the right to alter product specifications without prior notice. All rights reserved. © 2015 Don Whitley Scientific Limited.

Don Whitley Scientific Limited

14 Otley Road, Shipley, West Yorkshire, BD17 7SE, England. Telephone: +44 (0)1274 595728 Fax: +44 (0)1274 531197 Website: www.dwscientific.co.uk Email: sales@dwscientific.co.uk

Whitley Workstations H135 HEPA Hypoxystation

Your environment defined[®]

The H135 HEPA was developed for the study of cells under hypoxic and anoxic atmospheric conditions. With such accurate control and the ability to manipulate cells in situ without altering the incubation parameters, cell biology research can be performed over a comprehensive range of oxygen tensions with precision and total confidence.

The tallest, deepest, widest cell culture workstation in the current Hypoxystation range, the H135 HEPA has a huge internal capacity of 826 litres. The additional height and depth provide the ideal environment for the accommodation of a variety of items of laboratory equipment such as live cell imaging devices, microscopes and plate readers.

This Hypoxystation comes complete with internal HEPA filtration, environmental control systems, removable front, bespoke trolley and CO2 monitoring, full colour touchscreen interface and our unique oval glove port system that is fitted with sleeves for user comfort and ease of movement.

For users who wish to use a Seahorse XF Analyzer in their work, the H135 HEPA can connect to a Whitley i2 Instrument Workstation, which can house the Seahorse in an ideal, CO2 free working environment (as specified by the manufacturers). Plates can be passed from one chamber into the other via the Whitley Transfer Tunnel without subjecting cells to ambient conditions.

> For more information on the i2, please see our website or contact sales@dwscientific.co.uk



don whitley scientific

your environment defined[®]

O2 and CO2 levels monitored and adjusted continually to provide vou with better control and more confidence in your results

Whitley Internal HEPA Filtration Model H135 - exceeds the requirements of ISO 14644 Class 3

Designed specifically for researching cell biology under strictly maintained normoxic, hypoxic and anoxic conditions

Precise control designed to establish and maintain the environmental conditions you require throughout the entire incubation and working areas

Gas mixing achieved instantly via a unique, fully integrated control system - rapidly create your selected environmental conditions and minimise bench space requirements

Optional maintenance-free Automatic Humidification System provides a sterile source of additional humidity

Easy, accurate, 2-point oxygen sensor calibration

Fully automated O₂ calibration routine



Go online to see our workstations in operation:

www.hypoxystation.co.uk

Easily removable front for equipment transfer and thorough cleaning



15cm Single Sample Entry System - Ideal to introduce individual samples quickly.



Oval Sleeved Port



12 litre rapid cycle airlock



Equipment Enclosure

Contains all instrumentation to control oxygen, carbon dioxide, temperature and humidity.

Chamber

High quality clear acrylic structure provides optical clarity. The chamber is annealed (heat treated) at least twice during and after manufacture the annealing process relieves the stresses induced as a natural consequence of machining, forming and polishing acrylic. The chamber provides generous incubation and working areas.

Portholes

Unique oval portholes fitted with isolators providing the ability to use only one port at a time saving time and conserving gas

Optional Single Sample Entry System

Individual samples can be introduced in seconds into the workstation.

Oxygen Sensor and Calibration

Automatic oxygen sensor calibration without taking the sensor out of the chamber.

whitley H135 HEPA hypoxystation

Removable Front

Easily removable front for equipment transfer, the introduction of bulk quantities of samples and thorough cleaning.

Touch Screen

Full-colour 7" touch screen display featuring intuitive software to access a sophisticated microprocessorbased operating system.

HEPA Filtration System

The Whitley Internal HEPA Filtration System ensures that particulate is swiftly

removed

Intelligent Airlock Equipment and sample transfer is very convenient using the 12 litre airlock.

Whitley Internal HEPA Filtration System

All the atmosphere in the chamber passes through the filter hundreds of times an hour, which ensures the chamber environment is cleaned quickly. Placing the filter inside the workstation prevents the filter becoming saturated with moisture. Levels of atmosphere cleanliness inside the workstation exceed the requirements of ISO 14644 Class 3.

Options and Accessories

Ao6102 Vacuum Take-off

Provides a method of aspirating used liquid media from inside the workstation in safety and without the risk of spillage or contamination.

Ao6104 Humidification System

Adds sterile moisture and maintains a sterile environment.

A06106 Double Internal Power Socket

A06109 Data Logging

Allows the data recording of environmental parameters inside the workstation - oxygen, carbon dioxide, temperature and humidity levels and cabinet internal pressure. The collected data can be downloaded via a USB interface to a memory stick and transferred to a PC for further analysis.

Ao6108 Data Logging with Oxygen Profiling

Allows the user to pre-programme different oxygen levels. The user can determine how long the workstation atmosphere remains at a particular oxygen level or adjust to higher or lower gas concentrations.

	A02945	Small Cable Gland
	A06128	Large Cable Gland
	Ao6103 Ideal for	Single Sample Entry System (15cm) quickly introducing small quantities of samples (situated on left
	A07208	Single Sample Entry System (9cm)
	A07211 Provides examina	Spotlight directional lighting to assist in the early identification of sample gro tion.
	A07218	Enhanced Biological Containment
	L00040	Multi-Point Temperature Mapping
	L00050	Single Point Temperature Mapping

In the interests of a policy of continuous product improvement the company reserves the right to alter product specifications without prior notice.

The airlock cycle takes only 60 seconds.

hand side panel).

wth or other specimen

Remote monitoring and control over a secure internet link



Oxygen profiling allows the user to pre-programme different oxygen levels.



Convenient USB memory stick download of stored environmental parameters

